

870 KF Titrino plus



Water determination easy, safe and precise

870 KF Titrino plus – the KF Titrator with the best price-performance ratio

02

Intuitive and easy operation

The favorably priced 870 KF Titrino plus is Metrohm's Karl Fischer titrator for volumetric water determinations. It can be used to determine water contents from a few ppm up to 100% reliably and accurately in solid, liquid and gaseous samples. With its operating interface designed for routine users, the 870 KF Titrino plus is so easy to operate that it only requires a short familiarization period. Also its robustness makes it the ideal titrator for routine determinations.



The advantages of the 870 KF Titrino plus at a glance

03

- **Very simple installation**

The installation of the 870 KF Titrino plus is extremely simple. As soon as they are connected, the Exchange Unit, stirrer and compact USB printer are recognized and configured automatically.

- **Ready-to-use methods simplify your work**

You only have to decide whether you want to carry out a titer, blank or water determination and then start the titration by pressing just one key. Of course, you can also adapt the parameters yourself to your sample.

- **Formulas included**

The appropriate calculation formulas for each method are stored in the titrator. All you have to do is select the units for entering the sample weight and presenting the results; the 870 KF Titrino plus informs you about the necessary conversion factors.

- **Mouse or keyboard – choose your favorite!**

The 870 KF Titrino plus can be controlled either by mouse-click or keyboard – the choice is yours.

Highlights

- Straightforward installation
- Easy operation thanks to predefined methods
- Excellent precision thanks to the high-resolution measuring input
- Live curve for recognizing side-reactions
- Overflow protection for increased workplace safety
- Compact USB printer as an option





Precision you can rely on

Highest precision

Just like the Titrando high-end titrator, the KF titrator Titrino plus is also equipped with a high-resolution measuring input that guarantees results of the highest precision.

The live curve allows you to recognize side-reactions

Metrohm KF titrators are equipped with a display showing the live curve as standard. This allows you to recognize side-reactions in good time and take measures to suppress them.

Intelligent Exchange Unit for more comfort

The chip of the intelligent Metrohm Exchange Unit automatically provides the titrator with all the data (cylinder volume, type of KF reagent, validity of titer, etc.) that is necessary to carry out an error-free titration. This means that you always titrate under optimal conditions and, for example, do not have to repeat measurements because the determination has been carried out with an invalid titer.



Workplace safety is a priority

Overflow protection improve workplace safety

Great importance has been given to workplace safety: The new titrator allows safe handling of Karl Fischer reagents. The new safety parameter «Conditioning options» prevents the titration cell from overflowing during conditioning.

Environmentally friendly

Like all Metrohm KF titrators, the 870 KF Titrino plus can also be used with ethanol-based Karl Fischer reagents without any problems.

The 803 KF Titration Stand

The 803 KF Titration Stand is used for stirring and for exchanging the used working medium. The built-in pump aspirates off the used titration solution and adds new working medium without the cell having to be opened. This shortens the conditioning time and prevents the escape of solvents and their accumulation in the lab atmosphere.

05



870 KF Titrino plus – overview

Dosing elements	Exchange Units
Intelligent Exchange Units with integrated data chip	•
Steps per cylinder volumne	10'000
Operation, dialog	Keyboard and mouse
Stirrer, titration stand	803 KF Titration Stand 801 Magnetic Stirrer as an option
Stirrer connection	1 MSB connection (Metrohm Serial Bus, Daisy Chain)
Connections of balances, printer, USB keyboard, USB mouse, USB memory stick and barcode reader	Via 1 USB Slave Port, RS-232/USB Box (option)
Live curve display with graphic LCD	•
KF volumetric Karl Fischer titration with automatic conditioning	•
I _{pol} and U _{pol} – integrated programmable polarizer	•
Predefined methods including calculation formulars	•
Dialog language	Englisch, german, french, spanish, portugese, italian
Comprehensive GLP-/GMP funktions	•

Ordering information, options

- 2.870.0010

870 KF Titrino plus with one measuring interface, including comprehensive accessories with titration vessel and indicator electrode
- 2.803.0010

803 KF Titration Stand with stirrer and integrated membrane pump for reagent exchange
- 2.870.1010

870 KF Titrino plus with titration unit, indicaton electrode, 803 Titration Stand and 10 mL exchange unit
- 2.870.1020

870 KF Titrino plus with titration unit, indicaton electrode, 803 Titration Stand, 10 mL exchange unit and USB thermal printer Neo's
- 2.870.2010

Package including 870 KF Titrino plus, 860 KF Thermoprep and 803 Titration Stand

Intelligent Exchange Units equipped with data chip, including accessories and two buret tips, one of them with micro outlet valve.

- 6.3026.110

Buret volumen 1 mL
- 6.3026.150

Buret volumen 5 mL
- 6.3026.210

Buret volumen 10 mL
- 6.3026.220

Buret volumen 20 mL

Options for connecting a balance

- 6.2148.020

RS-232/USB Box
- 6.2151.030

Cable USB A – USB B

Optional compact USB printer

- 2.142.0100

Thermoprinter Custom Q3X



The optimal solution for every task



901 Titrande – volumetric high-end titrator for demanding KF titrations.

07



The 851 and 852 Titrande for water content determinations in the trace region. The picture shows the 851 Titrande.



The Karl Fischer oven technique for the easy preparation of difficult samples. The picture shows the 860 KF Thermoprep with manual sample preparation.



High sample throughputs are best dealt with by the fully automated 874 Oven Sample Processor. In the picture with 851 Titrande.

www.metrohm.com

