




PLANT MOISTURE SYSTEMS

Agriculture/Horticulture
Plant Growth Studies



Introduction

Skye have been manufacturing plant moisture systems since 1987. For nearly 20 years these instruments have been used to study the demand for water in plants by measuring the plant moisture potential or stress.

The demand for water in the plant is a combination of many factors including availability of water to the roots, environmental factors (for example temperature, wind, humidity, solar radiation) and the ability within the plant to move water.

The plant moisture potential measurement helps relate these factors to the state of the plant. It can be used to study the water need of the plant and the effects the local environmental conditions are having.

The plant moisture potential is measured using a Scholander type pressure chamber, otherwise known as a Plant Moisture System. These plant moisture systems have also been used for sap extraction studies. Skye Instruments Ltd has 2 different types of systems for this measurement and they are described below.



Analogue System

The low-cost analogue system, has a traditional dial gauge readout.

The water potential pressure of the plant sample is read off the large, clear dial gauge when moisture is first seen coming from the freshly cut stem.

The gauge has a maximum pressure needle which can be used to 'freeze' the reading, allowing the user to concentrate on the specimen and not on the gauge.

Digital System

The Digital System uses an accurate pressure transducer to give a digital readout via a clear display - in which the figures can be seen even in strong light.

Chamber Head

In collaboration with users of Plant Moisture System, we have 2 designs of head. The low pressure head has only one part and the seal is made by the gas pressure inside the vessel pushing against the soft rubber seal. This makes it very suitable for fleshy samples. The high pressure head, is recommended for more woody samples. Here the stem of the plant tissue is sealed in the neck of the vessel by a biconical compression seal.



Both heads are designed with a short seal length to accommodate 'awkward' samples with short petioles. Sealing pressure may be increased at will, even when the vessel is pressurised, by turning the seal clamp, should leakage occur.

The Skye system minimises damage to the stem by ensuring only the minimum pressure is applied to the petiole.

Gas Supply

In a laboratory the gas for pressurisation would normally be from large commercial cylinders of nitrogen. These would normally be rented together with a reduction valve from the local supplier. This is however not acceptable for portable use, and thus we offer for sale, small very portable compressed air cylinders with reduction valves. These may be refilled at thousands of centres worldwide, most commonly located in university towns and cities.





SPECIFICATIONS

	Analogue			Digital		
	SKPM 1405/40	SKPM 1405/50	SKPM 1405/80	SKPM 1400/40	SKPM 1400/50	SKPM 1400/80
Size 152 x 317 x 254mm	✓	✓	✓	✓	✓	✓
Chamber dimensions 135 x 70mm diameter 0.54 litre Volume	✓	✓	✓	✓	✓	✓
System weight - 9kg (Excludes gas bottle)	✓	✓	✓	✓	✓	✓
<u>Operating Pressure</u>						
0-40 bar	✓	✓	✓	✓	✓	✓
0-50 bar		✓	✓		✓	✓
0-80 bar			✓			✓
<u>Type of Chamber head suitable</u>						
Low Pressure (max pressure 40bar)	✓	✓	✓	✓	✓	✓
High Pressure		✓	✓		✓	✓
Dial Pressure Gauge (In units of Bar with max pressure reached indicator)	✓	✓	✓			
LCD Readout 12.7mm height, 2 decade ranges 0.01 Bar / 1KPascal resolution				✓	✓	✓

Supplied as complete systems with a range of sealing apertures
Accessories required : Gas Bottle and Regulator

SAFETY FEATURES

Skye's Plant Moisture System for the measurement of Plant Water Potential has a number of safety features included.

- The main chamber is machined from a single piece of solid, high grade brass. There are no metal joins and seams to fail under high pressures.
- The chamber head is secured using a quick fastening, multi start thread. This design ensures the chamber cannot be pressurised unless the head is secure.
- Each Skye plant moisture system is individually tested to a pressure far greater than the working pressure by an independent test house. A certificate of this test is supplied with every system.
- There is a safety release pressure valve fitted to each system, this is set to automatically evacuate the chamber if the maximum working pressure is exceeded at any time.
- Included in the cost of a Skye system is a pair of safety spectacles and a clear perspex safety shield. These ensure that no leaf or petiole debris can escape and be blown into the user's eyes.

PUBLICATIONS

R. Alonso, S. Elvira, V. Bermejo, B.S. Gimeno. Does drought stress protect *Quercus ilex* from ozone effects? Different response of subsp. *Ilex* and subsp. *Ballota*. 2005. Workshop: Critical levels of ozone: further applying and developing the flux-based concept Obergurgl, Austria, Nov 15-19 2005. www.uni-graz.at/pphwww/mitarb/mike/obergurgl/abstracts/poster_1.pdf

Pflanzenphysiologische Übungen. Humboldt University of Berlin, Institute of Biology, Plant Physiology. April 2005. www.biologie.hu-berlin.de/~plantphys/lehre/scripts/komplex2.pdf

A. Dattola, M.A. Lo Gullo, A. Motisi. Effetto del Portinnesto sull'Architettura Idraulica del nesto in Pianta Adulte della Nettarina 'Adriana'. 2003. IV Convegno Nazionale sulla Peschicoltura Meridionale. www.unipa.it/medpeach/atti/IV_Convegno/poster/Dattola_et_al_2.pdf

A. R. Reddy, D. Sundar, A. Gnanam. Photosynthetic flexibility in *Pedilanthus tithymaloides* Poit, a CAM plant. 2003. Journal of Plant Physiology. 160, 75 - 80 (2003)

B. Osborne, K. Black, G. Lanigan, M. Perks, G. Clabby. Survival on the exposed Limestone pavement in the Burren : photosynthesis and water relations of three co-occurring plant species. 2003. Biology and Environment : proceedings of the Royal Irish Academy. Vol 103B. No. 3, 125-137 (2003)

I. De Luis, J. J. Irigoyen, M. Sánchez-Díaz. Vapour pressure deficit reduces the beneficial effect of elevated CO₂ on growth of N₂-fixing alfalfa plants. 2002. Physiologia Plantarum. Vol 116, Issue 4, pp 497-502 H. Bauer. 2000. [Http://info.uibk.ac.at/projects/embolie/linkuk.html](http://info.uibk.ac.at/projects/embolie/linkuk.html)

M.P. Perks, S. Monaghan, C. O'Reilly, B.A. Osborne, D.T. Mitchell. Chlorophyll fluorescence characteristics, performance and survival of freshly lifted and cold stored Douglas fir seeds. 2001. Annals of Forest Science. Vol 58. (2001) 225-235

P. Giorio, G. Sorrentino, R. d'Andria. Stomatal behaviours, leaf water status and photosynthetic response in field-grown olive trees under water deficit. 1999. Environmental and Experimental Botany. Vol 42, Issue 2, 95-104.

G.E. Jackson, J. Irvine and J. Grace. Xylem Acoustic Emissions and Water Relations of *Calluna Vulgaris* L. at two Climatological Regions of Britain. Plant Ecology. 1999. 140, 3-14.

S.C. Clifford, S.K. Arnd, J.E. Corlett, S. Joshi, N. Sankhla, M. Popp, H.G. Jones. The role of solute accumulation, osmotic adjustment and changes in cell wall elasticity in drought tolerance in *Ziziphus mauritiana* (Lamk.). 1998. Journal of Experimental Botany. Vol 49. No. 323. 967-977.

A. Repellin, A.T. Pham Thi, A. Tashakorie, Y. Sahseh, C. Daniel, Y. Zuily-Fodil. Leaf membrane lipids and drought tolerance in young coconut palms (*Cocos nucifera* L.). 1997. European Journal of Agronomy. Vol 6, Issues 1-2, Pages 25-33

G.E. Jackson, J. Irvine and J. Grace. Xylem Cavitation in two mature Scots pine forests growing in a wet and a dry area of Britain. Plant, Cell and Environment. 1995. 18, 1411-1418. (Measurements on pine needles)

G.E. Jackson, J. Irvine and J. Grace. Xylem Cavitation in Scots pine and Sitka Spruce saplings during water stress. 1995. Tree Physiology. 15, 783-790. (Measurements on pine shoots)

G.E. Jackson, J. Irvine, J. Grace and A.M.M. Khalil. Absciscic Acid Concentrations and Fluxes in droughted Conifer saplings. 1995. Plant, Cell and Environment. 18, 13-22. (Xylem collection)

D. Le Thiec, M. Dixon, J. Garrec. Distribution and variations of potassium and calcium in different cross sections of *Picea abies* (L) Karst needles and *Fagus sylvatica* (L) Leaves exposed to ozone and mild water stress. 1995. Annales des Sciences Forestieres. Vol 52 (5). 411 - 422.

USERS

UK

- Forestry Commission Research Station
- Many UK Universities
- Silsoe College, University of Cranfield
- Horticulture Research International, East Malling Research
- Scottish Crop Research Institute
- Macaulay Land Research Institute
- Cambridge University Farm


Overseas

- Agricultural Research Institute, Cyprus
- Danish Institute of Agricultural Science, Denmark
- Debrecen University, Hungary
- University College Dublin, Ireland
- Trinity College Dublin, Ireland

- Central Mining Research Institute, India
- CNR ISAFoM, Italy
- Università Di Roma La Sapienza, Italy
- Kangwon Forest Research Institute, Korea
- Research Institute of Pomology & Floriculture, Poland
- Université Autonoma de Madrid, Spain
- Dpto. Biología de Organismos y Sistemas, Spain
- Fundación CEAM Centro de Estudios Ambientales del Mediterraneo, Spain
- USDA-ARS, USA

THOUSANDS OF CUSTOMERS IN OVER ONE HUNDRED COUNTRIES

Skye is a family run company and since 1983 has been exporting instruments to nearly every country in the world. We pride ourselves on customer care and our flexibility when it comes to providing the customer with what they need.



Skye Instruments Ltd
21, Ddole Enterprise Park
Llandrindod Wells
Powys
LD1 6DF
UK

Tel: +44 (0)1597 824811

Fax: +44 (0)1597 824812

Email: skyeemail@skyeinstruments.com

Web: www.skyeinstruments.com



**environmental and botanical
instrumentation**