

Lyssy L80-5000

Water Vapor Permeability Tester

Easy and reliable testing for water vapor permeability



The L80-5000 is the latest generation of the proven L80 series, which has been used successfully for decades in many labs.

With the L80-5000, the testing of water vapor permeability becomes very simple. The L80-5000 is easy to operate, and its reliable testing principle and low degree of maintenance have demonstrated to users worldwide the superiority of the Lyssy testing method.

The L80-5000 is easy to set up and operate through the built-in alphanumeric keyboard and digital display.

With a built in printer the L80-5000 also enables easy test documentation.

Features & Benefits

- Simplicity in operation due to the high degree of automation - the quality of tests performed is less operator dependent
- Easy to use operator interface
- Fast and accurate test result
- Extremely broad testing range, covering low and high permeability
- High repeatability of testing results
- Easy set-up of test parameters and sample data
- Complete traceability in test documentation, data logging and error reports
- Automatic equilibrium detection
- Automatic temperature control
- Easy to use test sample holders - no grease needed for sealing
- Built-in printer
- Low degree of maintenance


Sample cards

The preparation of a sample for the L80-5000 is accomplished in minutes. No grease or glue is required for a tight seal around the sample in the chamber. The test sample is affixed to the self-adhesive sample card, which is inserted into the L80-5000 test chamber, separating the upper and lower chamber. The Lyssy L80-5000 is able to measure high permeability materials using special sample reduction cards that decrease the surface area of the sample.

This reduction of the surface prevents the system from getting saturated, and therefore it becomes possible to dry down the upper chamber of the instrument and obtain a permeability measurement.



Example of measuring cycle L80-5000



Technical Specifications

General description

Dimensions	480H x 400W x 470D mm
Weight	Approximately 26kg
Measuring range	0.03 - 10,000 g/m ² /day as standard
Sensor life	More than 5 years under normal conditions
Voltage	230 VAC or 115/100 VAC
Conforms to the following standards	ASTM E398-03, ISO/CD 15106, JIS K 7129, TAPPI T523 om-82, NF H00-044

Measuring

Measuring temperature range	5-70°C. Practical range 30-70°C, with built-in temperature control. For measuring below 30°C an external cooling water thermostat is required
Humidity range	10% or 35% RH in measuring chamber (equivalent to 90% or 65% RH differential over the sample)

Sample Requirements

Measuring area	Low permeability samples - 50cm ² High permeability samples - 2.5cm ²
Sample thickness	Up to 12mm
Minimum sample size	10 x 10cm

Data logging

Memory	1 record with 1427 measurements or 109 records with 10 measurements or 353 records with 1 measurements
Printer	Built-in 80 mm thermal printer

User Interface

Keyboard	Alpha numeric
Display	Vacuum fluorescent display

Operational Environment

Ambient temperature	5-40°C
Ambient humidity	10-90% RH (non-condensing).

Fast, accurate and versatile

The L80-5000 is very fast at measuring, and the more permeable the sample the shorter the measuring time. In addition to its short testing times and broad testing range, the L80-5000 has a high degree of accuracy.

A series of tests by users have proven that the standard deviation of the Lyssy instruments is lower than $\pm 5\%$, and the reproducibility tolerance is as low as 1.5% (depending on the material).

The high accuracy and extremely broad testing range of 0.03-10,000 g/m²/day are achieved by using a very sensitive and reliable humidity sensor, which is located directly in the measuring chamber. This test method is the best reproduction of real-life conditions, since no carrier gas or extractive measuring technique is used.

The humidity sensor is very stable, regardless of the humidity range in which it is operated. As a result, the L80-5000 alternates easily and quickly between low and high permeability measurements - typically, the change can be done in one hour or less using the function called "Simulated Test standard". The tester automatically detects the attainment of equilibrium when the sample has stabilised.

Measuring temperature

The water vapor permeability of many products is strongly temperature dependent. That is why the L80-5000 can be used at precisely the required measuring temperature. However in order to test at temperatures below 30°C, the measuring chamber must be cooled to the required temperature via an external water cooling thermostat, which can be supplied as an optional accessory.

Illinois Instruments, Inc (U.S)

2401 Hiller Ridge Road
Johnsburg, Illinois 60051
U.S.A
Tel: +1 815 344 6212
Fax: +1 815 344 6332
E-mail: sales.usa@systechillinois.com
www.systechillinois.com

Systech Instruments Ltd (UK)

17 Thame Park Business Centre,
Wenman Road,
Thame, Oxfordshire OX9 3XA
Tel: +44 (0)1844 216838
Fax: +44 (0)1844 217220
E-mail: sales.uk@systechillinois.com
www.systechillinois.com

Illinois Instruments (Thailand)

6th fl Nopnarong Bldg No7
Ladprao23, Jatujak, Bangkok 10900
Thailand
Tel: +66 (0)2938 0798
Fax: +66 (0)2938 1058
E-mail: sales.ap@systechillinois.com
www.systechillinois.com

Systech Illinois (China)

Room 1105, Forte Building
No. 910 Quyang Rd, Hongkou district,
Shanghai, China 200434
Tel: +86 21 65533022
Fax: +86 21 65539651
Email: info@systechillinois.cn
www.systechillinois.cn