

New
in the
team?

YEP!



OMNIS MEETS CUBIS® II

SEAMLESS INTEGRATION OF THE LABORATORY BALANCE
FROM SARTORIUS



OMNIS IS NOW ALSO AN ANALYTICAL BALANCE

Thanks to an extension (QApp) implemented on the balance, CUBIS® II from Sartorius has become a fully integrated member of the OMNIS platform. Be it compliance with regulatory requirements, process design, or network integration – the integration of the high-end laboratory balance from Sartorius in OMNIS meets superior standards regarding performance and usability.



Data security and integrity

CUBIS® II not only sends the measurement results to OMNIS, but also the complete metadata for every weighing-in or weighing-out performed on the balance. When integrated in a local network, encryption of any data exchanged makes sure the highest standards of data integrity and security are met at all times.

Central user management

The user administration for the CUBIS® II balance is managed and verified centrally in OMNIS. The complete Audit Trail of all events on the balance is recorded in OMNIS. All requirements regarding regulatory compliance (FDA and EurdraLex) are fully met.

Maximum efficiency

The weighing-in or weighing-out of a sample or even of a complete sample series can be started right from the balance – without the need to even start OMNIS. Tolerance limits regarding sample weights are graphically represented on the display of the balance. When part of a local network, several OMNIS Clients can «share» the same CUBIS® II balance.

