

thermal hazard technology

Macro Reaction Calorimeter MRC™

The cost effective alternative



Base unit plus:

**Reactor vessels to 2 litres
Pressure to 10 Bar
Automated reagent
addition
Reflux and Distillation
Excellent control and
sensitivity**

**Inexpensive and versatile
Reaction Calorimetry for all!**

**1 North House, Bond Avenue, Bletchley
MK1 1SW, UK**

**Tel: +44 1908 646800 Fax: +44 1908 645209
Info@thermalhazardtechnology.com**

www.thermalhazardtechnology.com



Flexible and easy to use yet highly accurate. The THT MRC™ will at a fraction of the traditional price provide process reaction calorimetry for all.

The price conscious nature of this product will enable reaction calorimetry to be performed by process development groups as well as the more usual process safety groups. This will allow for faster development of process from bench to pilot and full plant.

Additional options for this system provide dosing, pressure handling and reflux reactions to be studied.



**1 North House, Bond Avenue, Bletchley
MK1 1SW, UK**

**Tel: +44 1908 646800 Fax: +44 1908 645209
Info@thermalhazardtechnology.com**