

Our new sample collector
FRL 6/7/854

with a smaller footprint
for optimized laboratory space



ERWEKA

ERWEKAs compact

FRL 6/7/854

The FRL 6/7/854 combines a space-saving footprint with all the familiar functions of the large FRL generation. The legs have been moved inwards, so that the shortened FRL can stand safely on the laboratory table. It can be tested with 6, 7 or 8 channels and, depending on the application, can be equipped with different racks for different test tubes or HPLC vials.

Highlights

- 6.5 cm shorter and 5 cm lower than predecessor series, which saves 25 cm when rotated by 90°
- Works with IPC and PVP pumps
- Reduced leg distance for more stability on laboratory tables
- Easy cleaning and access to sample tubes
- Up to 26 rows for extended release formulation testing
- Special procedure for foaming media handling
- Automated flush routine between sampling intervals
- Different racks available (interchangeable with FRL x24 series racks)

Technical Data

Channels	6, 7 or 8
Max. capacity (not for 25 ml test tubes)	26 sample intervals
Valve	Integrated 3 way valves
FRL Racks	26x8 / 1.5 ml HPLC vials 26x8 / 4.0 ml HPLC vials 26x8 / 10 ml test tubes 18x8 / 25 ml test tubes
Weight	25 kg
Dimensions: width/depth/height	700 mm/515 mm/585 mm
Voltage $\pm 10\%$	115 - 250 VAC; 50/60Hz



ERWEKA GmbH

Pittlerstr. 45
63225 Langen
Germany

E-Mail: sales@erweka.com
Phone: +49 6103 92426-200
Fax: +49 6103 92426-999

Technical specifications of products described are stated without warranty and subject to change at any time without further notice. **v.1.2.4.19**

ERWEKA