



TacticID

B&W TEK

Now a Metrohm Group Company

Handheld Raman Analyzers for Rapid Identification of Suspicious Materials


The TacticID family from B&W Tek is a complete group of field-ready handheld Raman analyzers designed for easy ID of the greatest variety of hazardous materials, explosives, chemicals, and narcotics.

- Able to measure dark and colored samples with optional accessories
- Non-destructive testing of solids, liquids, gels, and powders
- Critical component mixture analysis can identify threats in mixed samples



Applications

- Identification of bulk explosive materials and hazardous chemicals
- On-site identification of narcotics and controlled substances
- Non-destructive forensic evidence identification
- Identification of precursors, binders, and cutting agents



Features

Ease of Use: Designed for technical and non-technical users alike, featuring simple, easy to understand menus and commands, and operation using software or hardware buttons.

Onboard Camera: Capture evidence photos to accompany scan reports.

Rugged: Up to IP68 rated and U.S. MIL-STD-810G certified for drop testing.

Extensive Libraries: Over 10,000 on-board library items available - explosives*, chemicals, and narcotics with user creation and import options

Easy-to-Read Results: Color-coded Red, Yellow, and Green results with chemical safety symbols.

Bright Touch Screen: Ultra crisp HD display with touch screen capability for on-the-spot note and text inputs.

Designed for Safety: Scan delay, adjustable laser power, and hands-free operation options.



For More Information: www.tacticid.com